

code rings & tapen tell



CODE RINGS

For colour coding hand instruments.


- Silicon Rubber rings, fit all hand instruments.
- Sterilizable up to 240°C / 464°F.
- In 9 different colours for easy identification.

CODE RINGS ROULETTE 225 RINGS

9 colours x 25 rings

BOX OF 50 RINGS - ONE COLOUR.

Black, Blue, Brown, Green, Grey, Orange, White, Yellow and Red.



TAPEN TELL

- A surgical and dental colour coding tape in 8 colours.
- Spool of 3200 x 7 mm..
- Autoclavable.
- Black
Blue
Brown
Green
Orange
White
Yellow
Red

**FOR EASY MARKING AND IDENTIFICATION
OF HAND INSTRUMENTS**



— SWEDISH DENTAL SUPPLIES AB —

Södervägen 30, SE-232 52 Åkarp, Sweden

E-mail: info@swedent.se, Fax +46 40 46 05 13, Phone +46 40 46 02 01, www.swedent.se