

screw posts

gold plated
or titanium



Root canal screws conical shaped to avoid unwanted wedge effect. Electrolytically plated with a thick layer of 24k gold. Available in 6 sizes and 4 length.

Size	1	2	3	4	5	6	
S	●	●	●	●	●	●	Short 8,0 mm
M	●	●	●	●	●	●	Medium 9,5 mm
L	●	●	●	●	●	●	Long 12,0 mm
XL				●		●	Extra long 14,0 mm
∅	1,05	1,20	1,35	1,50	1,65	1,80	Dimension

● ASSORTMENT BOX OF 60 PCS
GOLD PLATED or TITANIUM

3 pcs of 20 sizes, cross key and hallow key

● ASSORTMENT BOX OF 120 PCS
GOLD PLATED or TITANIUM

6 pcs of 20 sizes, cross key and hallow key

● ASSORTMENT BOX OF 240 PCS
GOLD PLATED or TITANIUM

12 pcs of 20 sizes, cross key and hollow key

● REFILL BOX GOLD PLATED

Box of 12 pcs one size

● REFILL BOX GOLD TITANIUM

Box of 6 pcs one size



— SWEDISH DENTAL SUPPLIES AB —

Södervägen 30, SE-232 52 Åkarp, Sweden
info@swedent.se www.swedent.se